

## Aquaexcel repository - BioWes

### Installation

1. Connect to Aquaexcel server - <http://160.217.215.252/> (put it to your browser)
2. Click to Registration
  - a. Type your email, first name, last name and create password
  - b. Your email and password are your login information to the system
  - c. The request is sent to the repository administrator (Petr Cisar)
  - d. The administrator will create the account – you will get the confirmation by email
3. Install Protocol manager to local computer
  - a. Download installation of Protocol Manager - <http://download.datapartner.cz/data/BioWes/BioWes.msi>
  - b. Run the file
  - c. Install Protocol manager to the location where you have full control of the files (create, delete, modify)
  - d. Test the system using predefined account
    - i. The repository contains predefined testing account with predefined protocols and data for demonstration purposes
    - ii. For more details see the Testing manual

### Using Protocol manager

1. Run Protocol manager – BioWes.exe
2. Login using your credentials (obtained from new user creation)
3. Design protocol template
4. Realize real experiment
5. Download data from server back to computer

### Using web interface

1. Put <http://160.217.215.252/> to your browser
2. Login to the server using your credentials (obtained from new user creation)
3. See details about experiment
4. Share experiment
5. Search for particular experiment

**For more information see the video guides.**